

## Expression of proteins in *E.coli* with pBAT vectors

A general protocol for producing protein in *E.coli* using T7lac promoter based expression vectors (or any other vector which is controlled by *lac* repressor).

### Protocol

1. Start o/n culture of your construct in BL21(DE3) strain. Use always cells from a fresh plate, never from directly from glycerol stock. If possible plate the cells on an LB plate the day before. If you have an unstable plasmid, you might wish to skip the overnight liquid culture and start expression directly from the plate. Prepare 10 ml for each liter you plan to express the next day + some extra. Grow the cells in LB supplemented with ampicillin (100  $\mu\text{g}/\text{ml}$ ) at 37°C. If you have plasmid encoding the rare Arg-tRNA (pUBS-520) in the cells, add also 25  $\mu\text{g}/\text{ml}$  of kanamycin.
2. Put your LBs into 37°C to have them warm next morning.
3. In the morning pour 1 liter of LB or 2TY into 2 l shaker flask. Add the antibiotics as above and then 10 ml of the overnight culture. Or if you start directly from a plate, make a smaller starter culture with few colonies and add that to the final volume after an hour or two.
4. Grow at 37°C until OD<sub>600</sub> is 0.7 - 1.0 (usually 3 - 4.5 hours)
5. Add 1 ml of 400 mM isopropyl- $\beta$ -thiogalactopyranoside (IPTG) into each flask to induce the expression. If you wish to do the expression at lowered temperature (which can increase the solubility of your protein), transfer the flasks to desired temperature and allow sufficient time for cooling before induction.
6. Express the protein for 3 hours at 37°C. At lower temperatures the expression can be continued longer — 4 hours at 25°C and 20 hours at 15°C.
7. Spin the cells down 5000 rpm 10-15 minutes.
8. Carefully pour supernatant back to flask and add Virkon to decontaminate them. Let the bottles stand upside down to remove the last drops of media.
9. Add 20 ml of your lysis buffer to one bottle and vortex to resuspend the cells. Pour the suspension to next bottle and vortex again. Resuspend all the cells this way. In the end add few milliliters of fresh buffer to the first bottle, collect residual cell suspension by pipetting and do the same for the rest of the bottles. Combine all and either freeze at -20°C or proceed with lysis — French press or sonication.

## Solutions

### LB media

5 g Yeast extract  
10 g Bacto tryptone  
10 g NaCl  
H<sub>2</sub>O to 1 l, autoclave and store at room temperature

### 2YT media

10 g Yeast extract  
16 g Bacto tryptone  
5 g NaCl  
H<sub>2</sub>O to 1 l, autoclave and store at room temperature

### 100 mg/ml ampicillin stock (1000x)

1 g of ampicillin  
Sterile H<sub>2</sub>O to 10 ml  
Filter sterilise with 0.2  $\mu$ m filter  
Store in aliquotes at -20°C

### 25 mg/ml kanamycin stock (1000x)

250 mg kanamycin  
Sterile H<sub>2</sub>O to 10 ml  
Filter sterilise with 0.2  $\mu$ m filter  
Store in aliquotes at -20°C

### 400 mM IPTG stock (1000x)

953 mg IPTG  
Sterile H<sub>2</sub>O to 10 ml  
Filter sterilise with 0.2  $\mu$ m filter  
Store in aliquotes at -20°C